**SORIMEX**

**PL Elektrody EKG TAB**

**EN ECG TAB electrodes**

**NO EKG-elektroder TAB**

SORIMEX sp. z o.o.  
Równinna 25, 87-100 Toruń, Poland  
office@sorimex.eu ♦ www.sorimex.eu

Rev. 2026-04-24

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PL

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EN

## 1. Intended use

The electrode is applied to the surface of the body to transmit the electrical signal at the body surface to a recording and monitoring device that produces an electrocardiogram or vectorcardiogram.

## 2. Patient indication

Adult & pediatric.

## 3. Contraindications

There are no known contraindications for use of the device.

## 4. Application instructions

- 4.1. The electrode site should be clean, dry and free of skin oil before application of the electrode to allow optimal electrode adhesion to skin. Remove any excessive hair by shaving if required.
- 4.2. For best trace quality, skin abrasion should be part of the skin preparation. Use Abrasive Pads as part of the skin preparation.
- 4.3. Carefully peel off the electrodes one by one from the protective foil and place them on the designated monitoring locations on the skin. After removing all the electrodes, the protective foil should be immediately disposed of in the appropriately labeled waste bin.
- 4.4. Apply the electrodes with the connector facing downward and firmly press the electrodes to ensure a uniform contact with the skin.
- 4.5. Place all the electrodes on the patient before connecting the connectors.
- 4.6. Attach the metal part of the connector to the electrode in an area without gel. Ensure that the connector and the cable do not come into contact with the electrode gel.
- 4.7. Protect the remaining electrodes in the original bag to prevent moisture loss by securing the bag's opening.
- 4.8. Perform ECG procedure.
- 4.9. After conducting the examination, disconnect the connectors from the electrodes' cables.
- 4.10. Remove the used electrodes from the body and dispose of them properly.

## 5. Disposal

- Dispose off as hospital contaminated waste, in accordance with regulations in force.
- Recycle package.

## 6. Caution

- The electrodes should only be used by or in consultation with a health care provider familiar with their proper placement and use. Any allergic or other skin reactions should be immediately evaluated by a physician. For use only on intact, clean skin (do not use over open wounds, lesions, infected or inflamed areas).

- The product is not for internal use. In case of accidental ingestion of the electrode consult a physician immediately. The liner card should be immediately disposed off in the assigned waste bin for disposal after removal of all the electrodes.
- For single use only. Use on other patients can cause cross infection. Do not soak, rinse, or sterilize this device, as these procedures may leave harmful residues or cause malfunction of the device. The design and materials used are not compatible with conventional cleaning and sterilization procedures. Careless removal of electrode may cause damage to the skin. Should the electrode be difficult to remove, use water to dissolve the adhesive.
- US: Rx only

## 7. Precautions

- Self-adhesive electrodes should be replaced if they no longer stick firmly to the skin.
- During surgical procedures, electrodes should be placed as far as possible from unwanted RF current flow.
- Leaving the packaging open can cause the gel on the electrodes to dry out, so it is recommended not to open it unless necessary.
- Electrodes with no lead wire attached should only be used with lead wire in compliance with IEC 60601-1 and electrodes must only be used with ECG monitoring devices in compliance with IEC 60601-1.

MR Conditional electrodes only: It is important to be aware that the MR Conditional marking applies only to the electrodes and not to any monitors, cables or lead wires which may be attached to the electrode. Place the electrodes and electrode lead wires according to local safety guidelines for recording ECG in a MRI environment.

The MR Conditional electrode shall only be used with lead wires complaint with use in a MRI environment.

- Magnetic Resonance Conditional. Static magnetic field of 1.5 Tesla and 3 Tesla only. Maximum spatial gradient field of 25.000 Gauss/cm or 250 Tesla/m. Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning.

## 8. Potential adverse events

Potential adverse events in relation to the use of surface electrodes are skin reactions such as skin burn, skin redness, skin itching and small blisters.

## 9. Duration of use

Electrodes are applied in contact with skin surface for up to 30 minutes.

## 10. Storage

- Store in original packaging protected from weather and within the indicated temperatures.
- Once the pouch is opened, use electrodes promptly. Protect remaining electrodes in original pouch from moisture loss by folding-over the top and securing.
- Remaining electrodes should be used within 30 days of opening the pouch.

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

## 11. Explanation of symbols

CE mark (conformance to the Regulation (EU) 2017/745)

Product does not contain natural latex

Medical device

Product does not contain polyvinyl chloride

Do not re-use

Product has an Ag/AgCl sensor

Read the instruction for use

Caution

Temperature limit

Do not use if the package is damaged

Protect against sunlight

Manufacturer

Keep dry

MR conditional

Catalogue number

Indicates the date

Product index

Production batch number

Product UDI code

Pieces per packaging

## 12. Product incident reporting

Any serious incident related to the device should be reported to the manufacturer and to the authority competent for the place of residence of the user or patient.

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